Healthometer®

Professional Dial Scale

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Báscula Profesional

International Operating
Instructions

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Introduction

Congratulations! You have just purchased one of our finest professional scales. It is the latest in a line of precision weighing instruments that Health o meter has been building since 1919.

Ourexperttechnology, qualityand innovation haseamed us the confidence of doctors and the trust of consumers to bring medical excellence into your home.

This scale has been engineered to give you a lifetime of service. If you haveanyquestionsconcerning this scale, please refer to the following information:

If you require additional information, please call us at: 1-800-638-3722.

Operating the Scale

- 1 Make sure the scale is on a level surface. While the unit is designed with feet for weighing on carpet, a very heavy shag rug, or extremely plush carpet could affect performance.
 - 2 Step on and off the scale a few times before attempting to weigh yourself. This Will align the internal mechanism.
 - 3 Adjust the pointer to zero, if necessary, using the zero adjust knob at the top of the scale.
- Note: As you use the scale over a period of time you may flnd the pointer “drift” sllightly from zero. The best time to check the zero adjustment is after weighing when you step off the scale. If the pointer does not return exactly to zero, readjust with the zero adjust knob.**
- 4 Step on scale with your weight equally distributed on both feet while standing still.

Weighing Hints

- 1 Weigh yourself at approximately the same time each day. The foods you eat, liquids you drink and your level of activity Will cause your weight to vary at different times throughout the day.
- 2 For more realistic results, weigh yourself without clothing. The average weight of clothing, including shoes, is about 2-5 lbs. for women and 5-8 lbs. for men.

Maintaining the Scale

- 1 The baked enamel finish, lens and mat can be cleaned with a mild soap and damp, soft cloth. **DO NOT USE ABRASIVE OR STRONG CHEMICAL CLEANERS.**
- 2 All bearings and working parts are permanently lubricated. **DO NOT USE OIL OR LUBRICANTS.**
- 3 If your scale fails to work properly, do not attempt to dismantle it. Please refer to the Warranty Section of this booklet for servicing instructions.

Safety Precautions

- 1 Before using this scale, please read all instructions.
- 2 Health o meter scales are designed for weighing people, and are not appropriate for commercial or industrial use. This is designated by the words "Not Legal for Trade" on the scale dial.
- 3 This scale should not be used in the presence of flammable fumes.
- 4 This unit should not be used in lieu of professional diagnostic devices-consult your physician.

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LIFETIME LIMITED WARRANTY

This Health o meter® scale is warranted to be free of manufacturing defects for the Lifetime of the original purchaser or, if a gift, the original recipient. This warranty does not cover damage caused by accident or misuse.

If your scale becomes defective within this period, Health o meter agrees to repair or replace the product at its discretion free of charge and Will return the product postpaid.

In order to obtain warranty service, send your scale to Health o meter, attention Technical Services. Include a note giving the reason for returning your scale and print your name, address and telephone number. Make certain that your scale is properly packaged to avoid damage in transit to the factory. Please allow four weeks for service.

If the scale is not covered by warranty, or has been damaged, an estimate of repair costs or replacement cost Will be sent to you for approval prior to servicing or replacement.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state. Any implied warranty is limited in duration to the lifetime provided in this, the only, expressed warranty. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

In no event shall Health o meter be responsible for incidental or consequential damages, nor damage due to misuse or the use of any unauthorized attachment; nor damage caused by use of an electrical circuit not specified on the Health o meter product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

If you have any questions concerning the performance of this scale, please call Health o meter, Technical Services: 1-800-638-3722.

Technical Services Dept. Hours are 7:00 A.M. to 5:00 P.M. Central Standard Time.

WARRANTY APPLIES TO PRODUCTS SOLD IN U.S.A. AND CANADA OUTSIDE
U.S.A. AND CANADA, CONTACT DISTRIBUTOR FOR WARRANTY.

Heatthometer™

A Product of Signature Brands, Inc.
Bridgeview, Illinois 60455